


CranioStrip

Resorbable Inter-cranial Burrliner



KEY BENEFITS

The Osteopore CranioStrip is made of a bio-resorbable material with a inter-connected porous 3D architecture that allows for rapid and homogenous vascularisation after wound closure. This allows for new bone to fill up into the porous space ensuring the implant is well integrated with the surrounding calvarial bone.

1. FDA-cleared for craniofacial indications
2. Easy handling without need for external fixation
3. Excellent biocompatibility
4. Promotes bone formation

INDICATIONS

The Osteopore CranioStrip is indicated for use as a void filler between the cranial flap and the surrounding cranium after a craniotomy procedure.

CONTACT INFORMATION

Osteopore International Pte Ltd
10 Science Park Rd, The Alpha, #02-28, Singapore 117684
Email: skyew@osteopore.com.sg
Tel: +65 6779 4842
Fax: +65 6775 8625
www.osteopore.com.sg



Caution: Always refer to the package insert, product label, and/or user instructions before using any Osteopore product. The information presented in this brochure is intended only to demonstrate an Osteopore product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. This device can only be used upon prescription by a surgeon.

Osteopore

CS Nov 2008 Rev 1